# G3600/G3601 Inert Gas Oxygen Analyzing System

With Durable Zirconia Sensors



Perfecting Sensible Technology









## **Highest Quality at Lowest Maintenance**





G<sub>3600</sub> Inert Gas Oxygen Analyzing System—single board

#### The G<sub>3600</sub>/G<sub>3601</sub> Measure Oxygen in Inert Gas

The G<sub>3600</sub>/G<sub>3601</sub> Inert Gas Oxygen Analyzing System is a user-friendly and robust system for accurate and reliable oxygen measurements in safety critical applications.

#### **New Oxygen Analyzers are MED Approved**

The core of the system is the G<sub>36</sub> Oxygen Analyzer. The G<sub>36</sub> is approved under the European Marine Equipment Directive (MED), becoming the first system to be certified under the new heading A 1/3.54 for fixed oxygen analyzers.

Since the end of July 2009, fixed oxygen analyzers have to have MED type approval on European flagged ships – this includes most EU flags plus Norway.

#### **Obtain Redundancy with a Double Board**

The G<sub>3600</sub> is a single boarded system (see picture to the left) that fulfills all requirements and tasks at minimum cost.

The G<sub>3601</sub> is the double boarded version of the system (see picture to the upper right). This system offers redundancy for uninterrupted operation. If one sensor or analyzer fails, just switch to the other within seconds.

## Cost-Effective Redundancy with Double Board



### **Key Features**

- Certified under the M.E.D.
- Type Approved by DNV, BV, and Lloyd's Register











- Avoid delays
- Protect cargo, ship, and crew
- Easy to install, easy to maintain, easy to use
- Easy calibration by crew
- Option for simplified artificial calibration
- Configurable measuring range and signal outputs
- Configurable alarms
- Long time sensor stability easy replacement
- Inexpensive spare parts
- Graphic display interface via touch screen
- ✓ Worldwide customer support via service partners

#### **Durable and User-Friendly Zirconia Sensor**

The G<sub>36</sub> Oxygen Analyzers provide accurate, real-time measurements based on a new type of zirconia cell suitable for the harsh and stressful marine environment. The zirconia sensors are cost-effective and reliable. They are easily stored and have unlimited shelf-life.

#### **Easy Configuration**

The G<sub>36</sub> Oxygen Analyzers are easily set up to work seamlessly together with the inert gas generator. The interface allows you to freely configure both the analog output signal as well as the alarm levels via the touch screen.

#### **Optional Sample Flow Control**

The G<sub>36</sub> Oxygen Analyzers allow the integration of additional sensors for customizing the G<sub>3600</sub> Oxygen Analyzing System to your needs. For example, it is possible to add a flow sensor to monitor the flow of the inert gas sample.

#### **Panel Mounted Oxygen Analyzer**

Both the G<sub>3600</sub> single board and the G<sub>3601</sub> double board can be delivered with the G<sub>36p</sub> Oxygen Analyzer that is designed for panel mounting. This is especially useful where an existing panel-mounted oxygen analyzer is replaced.

## Specifications - G3600/G3601 IGOA System

## G<sub>3600</sub>/G<sub>3601</sub> Board Specifications

G<sub>3600</sub> single board  $H\times W\times D$ :  $600\times 500\times 134$  mm - Weight: approx. 12 kg without water and packaging  $G_{3601}$  double board  $H\times W\times D$ :  $610\times 790\times 134$  mm - Weight: approx. 20 kg without water and packaging

Sample quality Pressure: 0.05 to 1 bar - Flow: 2 to 8 l/min - Temperature: 0°C to 70°C

Sample manifold 3 ports - 1/8" BSP connection

Zero test gas  $\setminus$  E.g. 2% O<sub>2</sub> in N<sub>2</sub> supplied by owner - max. 10 bar - Reduction station - 1/8" BSP connection

Span test gas \ Instrument air according to ISO 8573-1.3.3.2 - max. 10 bar - Reduction station incl. filter - 1/8" BSP connection

Filter retention \ 95% of 1 \mu particles

### **G<sub>36</sub> Oxygen Analyzer + Sensor Specifications**

Certificates & approvals MED by DNV - DNV Type Approval - Lloyd's Register Type Approval - Bureau Veritas Type Approval - 🕻 🤅

Sensor technology Heated zirconia type sensor

Measurement range 0.0 ... 21.0%

Repeatability +/- 0.1 % of the measurement range Accuracy +/- 0.5 % of the measurement range

Response time 90% of F.S. in less than 10 sec. at sample flow 2 l/min

Power supply G<sub>36a</sub>: 100...230 VAC / 50...60 Hz.

Output signal 2 × 4...20 mA - range selectable. Default: A-out1: 0.0...25.0 % 02 / A-out2 not in use

Max. load  $600 \Omega / 24 \text{ VDC}$ 

Alarm functions 02 low or high - 02 high-high - System Fail
Alarm relays 4 relays, volt free, 5A 24 VAC/VAC

Interface Touch screen  $71 \times 39$  mm with trend graph display

Alarm log History and alarm logs on SD cards

Ambient temperature 0°C to 55°C

Location Safe area e.g. engine room
Analyzer casing Aluminum casing IP67









### **Optional Equipment**

Extension board Analyzer as specified above

for additional analyzer Dimension:  $H\times W\times D$ :  $610\times 290\times 130$  mm - Weight: Approx, 7 kg without packaging

Sample flow control Additional flow transmitter mounted on the board and integrated into the G<sub>36</sub>

Remote digital display 22...250 VAC/DC with 2 configurable alarm relays - Ambient temp.:  $-20^{\circ}$ C to  $60^{\circ}$ C - Panel cut-out:  $44.5 \times 91.5$  mm

G<sub>36p</sub> instead of G<sub>36a</sub> Power: 24 VDC - Ambient Temperature 0°C to 70°C - Enclosure IP55 if panel mounted

Panel cut:  $154 \times 73$  mm (W×H) - Front:  $178 \times 95$  mm (W×H) - Depth: 71 mm + cables

other specifications as above

Other optional equipment Pre-filter for sample gas, signal amplifier, and signal amplifier for logarithmic output

Specifications subject to change without notice



